

# VELCRO WRAPS

## OPTIONS AND MEASURING

### FARROW WRAPS

JOBST FarrowWrap garments are a range of adjustable compression devices that are available in a variety of styles, fabrics, sizes and compression classes.

Made with durable material, they are designed to help reduce limb swelling and maintain volume reductions, as well as promote healing of venous leg ulcers. They also provide the same strong support as multi-layer bandaging but in a simplified, intuitive design.

*Available in beige only*



	XS	SM	MED	LG	XL
C - calf circumference (widest point)	36-43cm	42-50cm	48-58cm	53-63cm	58-68cm
B - Ankle Circumference (widest point)	21-25cm	25-30cm	30-36cm	36-42cm	42-50cm
A-D Regular Length	28-35cm	30-37cm	31-39cm	33-41cm	33-41cm
A-D - Long Length	33-40cm	35-42cm	36-44cm	38-46cm	38-46cm

#### CIRCUMFERENCES


Left Right

B \_\_\_\_\_

C \_\_\_\_\_

#### LENGTHS

A-D \_\_\_\_\_



### COMPREFLEX WRAPS

Lightweight, breathable, and easy to slide on the SIGVARIS CompreFLEX Transition (formerly Compreflex Lite) works well for patients who need higher compression and containment but need a product that is easier to get on than compression stockings. It also has a thicker liner with compression in the foot area.

*Available in black or beige*



#### CIRCUMFERENCES LEFT RIGHT

D \_\_\_\_\_

C \_\_\_\_\_

B1 \_\_\_\_\_

B \_\_\_\_\_

L \_\_\_\_\_

#### LENGTHS taken on lateral aspect

FIBULAR HEAD

L1

BK LENGTH

REGULAR: 30-33.9cm

TALL: 34-42cm

LATERAL MALLEOLUS (outer ankle bone)

	SM	MED	LG	XL	XXL
D	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
C	29-39cm	34-44cm	40-50cm	46-56cm	54-64cm
B1	24-34cm	29-39cm	34-44cm	39-49cm	44-55cm
B	19-26cm	21-30cm	26-36cm	31-41cm	36-46cm

REGULAR: 30-33.9CM / TALL: 34-42CM

**COLOR**  
Black Beige  
(please circle)

Patient Name: \_\_\_\_\_ DOB: \_\_\_\_/\_\_\_\_/\_\_\_\_

**FAX ABOVE MEASUREMENTS WITH FACE SHEET TO:  
1-866-388-5430**